

MEDIVATORS®
SCOPE BUDDY®
Endoscope Flushing Aid



MEDIVATORS®
SCOPE BUDDY®
Endoscope Flushing Aid

Visual display
of endoscope
flushing time

User-selectable
time setting

One-touch
start/stop
operation

Simple and quick
tubing connections



Consistent, hands-free flushing of endoscopes

Reduces the risk of repetitive motion injuries

- Replaces repetitive and labor-intensive syringing of channels
- Hands-free operation reduces muscle fatigue
- Simple push button operation
- Provides constant fluid delivery through all endoscope channels

Shortens manual cleaning time

- Thoroughly flushes most endoscopes in 35 seconds
- Flushes multiple endoscope channels simultaneously
- Quickly and quietly delivers detergent and rinse solutions through endoscope channels

Enhances compliance with endoscope cleaning protocols

- Continuous, timed fluid delivery ensures consistent, repeatable endoscope flushing
- Audible and visual alerts upon completion
- Compatible with all flexible endoscopes including ultrasound and those with elevator wire channels
- Works with your existing sink; compact size accommodates countertop placement or wall mounting
- Improves manual cleaning throughput





Scope Buddy Dimensions

9.5 H x 6 W x 8 D (inches)
 (24.1 H x 15.2 W x 20.3 D (cm))

Unit Weight

10 lbs (4.5 kg)

Electrical Requirements

110 – 230 V, 50/60 Hz

Electrical Safety Certifications

UL Standard 61010-1
 CANADA/USA C22.2 No. 61010-1
 ETL and cETL listed

Fluid Delivery Rate

1 liter/minute

Flush Time

0 – 99 minutes

Timer Control

Digital display with touchpad controls

System Pressure

Built-in safety bypass

Cycle Completion Indicators

Audible and visual alerts

Unit Placement

Counter-top or wall-mounting
 (bracket included)

Warranty

One year

Scope Buddy Ordering Information

Model number	Description
EFA-US-G	North America
EFA-IN-G	International

Accessory Ordering Information

Part number	Description
78399-669	Complete Endoscope Connection Kit
75227-679	Graduated Cylinder

Detergent Ordering Information

Part number	Description
ML02-0106	Intercept Detergent (4 one-gallon bottles/case)



Combine Scope Buddy with Intercept® detergent for...

- Effective endoscope bacterial biofilm penetration and removal
- Superior removal of biological and organic soils
- Assured endoscope compatibility

MEDIVATORS, SCOPE BUDDY and INTERCEPT are registered trademarks of MEDIVATORS Inc.



MEDIVATORS Inc. Headquarters
 14605 28th Avenue North
 Minneapolis, MN
 55447-4822 USA
 Tel: +1.800.444.4729
 Fax: +1.763.553.3387

MEDIVATORS BV
 Sourethweg 11
 6422 PC Heerlen
 The Netherlands
 Tel: +31.45.5.471.471
 Fax: +31.45.5.429.695

MEDIVATORS Asia/Pacific Pte. Ltd.
 1 A International Business Park
 #05-01 Singapore 609933
 Tel: +65.6227.9698
 Fax: +65.6225.6848

MEDIVATORS Beijing Representative Office
 Room 708, 7th Floor Kaiheng Center, Block B
 No. 2 Chaoyangmennei Street
 Dongcheng District, Beijing China 100010
 Tel: +8610.6567.8446
 Fax: +8610.6567.8445

www.medivators.com

© 2012 MEDIVATORS Inc. P/N 50098-031A